EN Instruction Manual

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www.rossmax.com

Warranty Registration (Must be completed within 10 days of purchase)

Customer Name: _____ **City:** _____
Street Address: _____ **ZIP code:** _____
State: _____ **ZIP code:** _____
Telephone: _____ **E-mail address:** _____

Gender: Male Female **Age:** _____

Product Information

Date of purchase: _____

Store where purchased: _____

Price Paid (excl. Tax): _____

This instrument is covered by a 3 year guarantee from the purchase date. The guarantee is valid only on presentation of the guarantee card completed by the dealer conforming purchase date or the receipt. Nebulizer Components are not included. Opening or altering the instrument invalidates the guarantee. The guarantee does not cover damage, accident or non-compliance with the instruction manual. Please contact Rossmax Service or distributors.

ISO 9001

ISO 14001

ISO 13485

ISO 27001

ISO 22700

ISO 22703

ISO 22704

ISO 22705

ISO 22706

ISO 22707

ISO 22708

ISO 22709

ISO 22710

ISO 22711

ISO 22712

ISO 22713

ISO 22714

ISO 22715

ISO 22716

ISO 22717

ISO 22718

ISO 22719

ISO 22720

ISO 22721

ISO 22722

ISO 22723

ISO 22724

ISO 22725

ISO 22726

ISO 22727

ISO 22728

ISO 22729

ISO 22730

ISO 22731

ISO 22732

ISO 22733

ISO 22734

ISO 22735

ISO 22736

ISO 22737

ISO 22738

ISO 22739

ISO 22740

ISO 22741

ISO 22742

ISO 22743

ISO 22744

ISO 22745

ISO 22746

ISO 22747

ISO 22748

ISO 22749

ISO 22750

ISO 22751

ISO 22752

ISO 22753

ISO 22754

ISO 22755

ISO 22756

ISO 22757

ISO 22758

ISO 22759

ISO 22760

ISO 22761

ISO 22762

ISO 22763

ISO 22764

ISO 22765

ISO 22766

ISO 22767

ISO 22768

ISO 22769

ISO 22770

ISO 22771

ISO 22772

ISO 22773

ISO 22774

ISO 22775

ISO 22776

ISO 22777

ISO 22778

ISO 22779

ISO 22780

ISO 22781

ISO 22782

ISO 22783

ISO 22784

ISO 22785

ISO 22786

ISO 22787

ISO 22788

ISO 22789

ISO 22790

ISO 22791

ISO 22792

ISO 22793

ISO 22794

ISO 22795

ISO 22796

ISO 22797

ISO 22798

ISO 22799

ISO 22800

ISO 22801

ISO 22802

ISO 22803

ISO 22804

ISO 22805

ISO 22806

ISO 22807

ISO 22808

ISO 22809

ISO 22810

ISO 22811

ISO 22812

ISO 22813

ISO 22814

ISO 22815

ISO 22816

ISO 22817

ISO 22818

ISO 22819

ISO 22820

ISO 22821

ISO 22822

ISO 22823

ISO 22824

ISO 22825

ISO 22826

ISO 22827

ISO 22828

ISO 22829

ISO 22830

ISO 22831

ISO 22832

ISO 22833

ISO 22834

ISO 22835

ISO 22836

ISO 22837

ISO 22838

ISO 22839

ISO 22840

ISO 22841

ISO 22842

ISO 22843

ISO 22844

ISO 22845

ISO 22846

ISO 22847

ISO 22848

ISO 22849

ISO 22850

ISO 22851

ISO 22852

ISO 22853

ISO 22854

ISO 22855

ISO 22856

ISO 22857

ISO 22858

ISO 22859

ISO 22860

ISO 22861

ISO 22862

ISO 22863

ISO 22864

ISO 22865

ISO 22866

ISO 22867

ISO 22868

ISO 22869

ISO 22870

ISO 22871

ISO 22872

ISO 22873

ISO 22874

ISO 22875

ISO 22876

ISO 22877

ISO 22878

ISO 22879

ISO 22880

ISO 22881

ISO 22882

ISO 22883

ISO 22884

ISO 22885

ISO 22886

ISO 22887

ISO 22888

ISO 22889

ISO 22890

ISO 22891

ISO 22892

ISO 22893

ISO 22894

ISO 22895

ISO 22896

ISO 22897

ISO 22898

ISO 22899

ISO 22900

ISO 22901

ISO 22902

ISO 22903

ISO 22904

ISO 22905

ISO 22906

ISO 22907

ISO 22908

ISO 22909

ISO 22910

ISO 22911

ISO 22912

ISO 22913

ISO 22914

ISO 22915

ISO 22916

ISO 22917

ISO 22918

ISO 22919

ISO 22920

ISO 22921

ISO 22922

ISO 22923

ISO 22924

ISO 22925

ISO 22926

ISO 22927

ISO 22928

ISO 22929

ISO 22930

ISO 22931

ISO 22932

ISO 22933

ISO 22934

ISO 22935

ISO 22936

ISO 22937

ISO 22938

ISO 22939

ISO 22940

ISO 22941

ISO 22942

ISO 22943

ISO 22944

ISO 22945

ISO 22946

ISO 22947

ISO 22948

ISO 22949

ISO 22950

ISO 22951

ISO 22952

ISO 22953

ISO 22954

ISO 22955

ISO 22956

ISO 22957

ISO 22958

ISO 22959

ISO 22960

ISO 22961

ISO 22962

ISO 22963

ISO 22964

ISO 22965

ISO 22966

ISO 22967

ISO 22968

ISO 22969

ISO 22970

ISO 22971

ISO 22972

ISO 22973

ISO 22974

ISO 22975

ISO 22976

ISO 22977

Introduction
Thank you for purchasing a NE100 Compressor Nebulizer. With proper care and use, your nebulizer will provide you with many years of reliable treatments. This unit operates on standard AC power. Treatments are delivered quickly, safely and conveniently, making this unit ideal for all ages. We encourage you to thoroughly read this guide-book to learn about the features of your nebulizer. Your compressor nebulizer should be used under the supervision of a licensed physician and/or a respiratory therapist. Together with your physician and/or pharmacist, you can feel comfortable and confident knowing that you are obtaining the most effective inhalation treatments for your respiratory condition.

NOTE: Your nebulizer is intended for use in treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy. Please consult with your physician and/or pharmacist to determine if your prescription medication is approved for use with this nebulizer. For type, dose, and regime of medication follow the instructions of your doctor or respiratory therapist.

This device fulfills the provision of the EC directive 93/42/EEC (Medical Device Directive) and the European Standard EN 13544-1:2007+A1:2009 (Respiratory Therapy Device Unit - Part 1: Nebulizing systems and their components).

PLEASE read this manual carefully before use and be sure to keep this manual.

Cautions
Please use general safety precautions when operating your nebulizer. This unit should be used only for its intended purpose as described in this guidebook and with medications only under the supervision and instruction of your physician. Do not use the device in anesthetic or ventilator breathing circuits.

Product cautions

READ THE FOLLOWING BEFORE USING

- Do not use this unit while bathing or showering.
- Do not immerse the unit in liquid.
- Do not handle the unit of power cord with wet hands.
- Do not use the unit in liquid.
- Do not use while bathing.
- Do not reach for anything that has fallen into water immediately upon the unit.
- Do not use the unit if it has any damaged parts (including plug). If it has been submerged in water or dropped, promptly send the unit for examination and repair.
- The unit should not be used where flammable gas, oxygen or aerosol spray products are being used.
- Keep the air vents open. Do not place the unit on a soft surface where the openings can be blocked.
- If the medication cup is empty, do not attempt to operate the unit.
- If any abnormality occurs, discontinue to use until the unit has been examined and repaired.
- The unit should not be left unattended while plugged in.
- Do not tilt or shake the unit when in operation.
- Disconnect the unit from the electrical outlet before cleaning, filling and after each use.
- Do not use attachments unless recommended by the manufacturer.
- Do not disassemble or attempt to repair the unit.
- Do not use the device in anesthetic or ventilator breathing circuits.

Operating cautions

- Close adult supervision is highly recommended when the unit is used by children or invalids.
- Keep your eyes away from the output of medication mist.
- The maximum capacity of the medication cup is 5 ml and should not be overfilled.
- Do not use this unit while operating a vehicle.
- If any discomfort or abnormality occurs, stop using the unit immediately.
- Do not use the device if the air tube is bent.
- Pentamidine is not an approved medication for use with this device.

Storage cautions

- Do not store the unit in direct sunlight, high temperature or humidity.
- Keep the unit out of reach of small children.
- Always keep the unit unplugged while not in use

Cleaning cautions

- Check air filter, nebulizer, mouthpiece and any other optional component before each use. Dirty or worn parts should be replaced.
- Do not immerse the unit in water. It may damage the unit.
- Disconnect the unit from the electrical outlet before cleaning.
- Clean all necessary parts after each use as instructed in this guidebook.
- Always dispose of any remaining medication in the medication cup after each use. Use fresh medication in the medication cup.
- Do not store the air tube with moisture or medication remaining in the air tube. This could result in infection as a result of bacteria.

MEDICAL DISCLAIMER:

This manual and product are not meant to be a substitute for advice provided by your doctor or other medical professional. Don't use the information contained herein or this product for diagnosing or treating a health problem or prescribing any medication. If you have or suspect that you have a medical problem, promptly consult your doctor.

Durable periods are as follows, provided the product is used to nebulize 2ml of medication 2 times a day for 8 minutes each time at room temperature (25°C). Durable period may vary depending on usage environment.

Main unit	5 years
Nebulizer Kit	1 year
Air Tube, Mouthpiece	1 year
Air Filter	60 days
Adult and child masks	5 years

Product specifications

Power	AC 230V/50Hz or AC 220V/60Hz or AC 110V/60Hz
Power Consumption	≤ 130W
Sound Level	≤ 60 dBA (1 meter away from NE100)
Compressor Pressure Range	≥ 29 psi (200 kPa)
Operating Pressure Range	≥ 15 psi (103 kPa)
Operating Flow Range	≥ 3.5 lpm
Operating Temperature Range	10°C to 40°C (50°F to 104°F)
Operating Humidity Range	10 – 90% RH
Operating Atmospheric Pressure Range	700-1060 hPa
Storage Temperature Range	-20°C to 60°C (-4°F to 140°F)
Storage Humidity Range	10 – 90% RH
Dimension (L x W x H)	204.5mmx148mmx108mm (8.05x5.83x4.25 inches)
Weight	1360g (without accessories)
Medication Capacity	5ml(cc)
Particle Size (MMAD)	≤ 3.0µm
Average Nebulization Rate	≥ 0.3 ml/min
Standard Accessories	Nebulizer Kit, Air Tube, Mouthpiece, Filters (5pcs), Adult and child masks

^aSubject to technical modification without prior notice.
^bPerformance may vary with drugs such as suspensions or high viscosity. See drug supplier's data sheet for further details.

A. Product identification			
1. Nebulizer Kit	2. Nozzle	3. Angled Mouthpiece	
4. Air Filter	5. Child Mask	6. Adult Mask	
7. Air Tube	8. Air Tube	9. Air Filter Slot	
10. Air Intake	11. Power Cord	12. Power Switch	

B. Assembling your nebulizer kit
Follow the cleaning instructions in this guidebook under "Cleaning technique" prior to using your nebulizer for the first time or after it has been stored for an extended period of time.

REMEMBER: Always unplug the compressor and make sure the power-switch is turned to the "OFF" position before cleaning, assembling and before or after each use.

- Place the compressor on a flat, stable surface within reach.
- Gently twist and pull straight up on the lid of the nebulizer to separate into two parts (medication cup and cover).
- The stem inside the medication cup inserts into the tube of the nozzle.
- Add the prescribed amount of medication in the medication cup.
- Reassemble the nebulizer by carefully twisting the medication cup and cover together. Make sure that the two parts fit securely.

WARNING: The symbol on this product means that it's an electronic product and following the European directive 2012/19/EU the electronic products have to be disposed on your local recycling centre for safe treatment.

C. Operating your nebulizer

The nebulizer is operate at up to a 45° angle. If the angle is greater than 45°, no aerosol will be generated.

- Attach one end of the air tube connector to the air output.
- Carefully attach the opposite end of the air tubing connector to the stem at the base of the nebulizer kit.
- Attach the angled mouthpiece or mask to the top of the nebulizer.
- The capacity of the medication cup is 2-5 ml.

NOTE: A 30-minute interval is recommended after each use. The compressor will automatically shut off if it becomes overheated. If it when this happens, immediately:

- Press the power-switch to the "OFF" position.
 - Unplug the power cord from the outlet.
 - Allow the motor to cool for 30 minutes before using the nebulizer again.
- Before restarting the unit, make sure that the air vents are not obstructed.

Operating your nebulizer

After every use:

- Unplug the unit from the power source.
- Allow the unit to completely cool.
- Carefully detach the air tubing from the nebulizer and pour out any remaining medication.

4. Follow the cleaning procedures provided in this guidebook.

D. Cleaning procedures

Rinsing technique (performed after each treatment or before first use):

- Make sure that the power-switch has been turned to the "OFF" position and the unit has been disconnected from the power source.
- Disconnect the air tube from the nebulizer device.
- Gently twist and pull up the cover of the nebulizer kit to open and separate.
- Rinse the nebulizer kit and components with hot tap water.
- Dry with clean towels or completely air dry.
- Reassemble the nebulizer kit.

NOTE: For the first time cleaning or after the unit has been stored for an extended period of time, thoroughly clean all components, except the air tube. The nebulizer kit is dishwasher safe.

Cleaning the compressor

Wipe the compressor daily using a soft cloth.

NOTE: Any other form of cleaning or cleaning agents may damage the finish of the unit.

E. Changing the air filter

It is important to change the air filter approximately when the air filter turns gray. It is recommended to change air filter every 2 months.

- Remove the air filter cover by gently pulling forward.
 - Discard the gray filter.
 - Replace with a new, clean air filter.
 - Securely re-attach the air filter cover to the unit.
- NOTE:** Air filters cannot be cleaned or washed. Only NE100 air filters can be used. Do not substitute alternate material such as cotton. Do not operate without an air filter.

Troubleshooting

- If any abnormality occurs during use, please check and correct the following:
- Unit does not operate when power switch is pressed. Check the AC connection to the outlet.
 - No misting or low rate of misting:
 - Check that there is medication in the nebulizer cup.
 - Check the main unit if there is any physical damage.
 - Check the position of the nozzle inside the nebulizer.
 - Make sure that air tube and other components are properly attached.
 - Check the air filter and replace if necessary.

Protection against electric shock:

- Class II equipment

Type BF applied parts:


- III and II parts and accessories

Protection against harmful ingress of water and particulate matter:

- IP21


Degree of safety in the presence of flammable anesthetics or oxygen:

- No APAPG (not suitable for use in the presence of flammable anesthetics or oxygen).

Guidance and manufacturer's declaration-electromagnetic emissions			
The NE100 is intended for use in the electromagnetic environment specified below. The customer or the user of the NE100 should assure that it is used in such an environment.			
Immunity test	EC 60801 test level	Completion level	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	Class B	The NE100 uses RF energy only for its internal function. Therefore, as RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Class A	The NE100 is suitable for use in all establishments. Harmonic emissions IEC 61000-3-2
RF emissions CISPR 11	Class B	Class A	The NE100 is suitable for use in all establishments. Harmonic emissions IEC 61000-3-2
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliance		those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Guidance and manufacturer's declaration-electromagnetic immunity			
The NE100 is intended for use in the electromagnetic environment specified below. The customer or the user of the NE100 should assure that it is used in such an environment.			
Immunity test	EC 60801 test level	Completion level	Electromagnetic environment-guidance
Electrostatic discharge IEC61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst supply lines IEC 61000-4-4	±2kV for power supply lines ±2kV for input/output lines	±2kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV line-to-earth	±1kV differential mode Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips IEC 61000-4-11	<5% UTI>95% dip in short interruption 40% <UTI>80% dip in long variations on power supply	<5% UTI>95% dip in short interruption 40% <UTI>80% dip in long variations on power supply	Mains power quality should be that of a typical commercial or hospital environment. The use of the NE100 requires continued operation during power mains interruptions. It is recommended that the NE100 be powered by an uninterruptible power supply (UPS) or a battery.
Power IEC 61000-4-7	3AIm	3AIm	The NE100 power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UTI is the a.c. mains voltage prior to application of the test level.			
Guidance and manufacturer's declaration-electromagnetic immunity			
The NE100 is intended for use in the electromagnetic environment specified below. The customer or the user of the NE100 should assure that it is used in such an environment.			
Immunity test	EC 60801 test level	Completion level	Electromagnetic environment-guidance
Portable and mobile RF communications equipment should be used no closer to any part of the NE100 including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:			
Where P is the maximum output power rating of the transmitter in (W) according to the transmitter manufacturer and d is the distance in metres (m).			
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur if the NE100 is not in normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the NE100.			
b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 µV/m.			
NOTE1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) tele-phones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the NE100 is used exceeds the applicable RF compliance level above, the NE100 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the NE100.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 µV/m.			

Recommended separation distance between portable and mobile RF communications equipment and the NE100			
Conducted RF IEC 61000-1-4	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz
Radiated RF IEC 61000-1-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz
Recommended separation distance according to frequency of transmitter (m)			
Power (W)	d=1.2√P	d=1.2√P	d=1.2√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter. Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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