



Infrared Ear Thermometer

Model: RA600



- 1 second measurement or continuous scanning mode
- Ear temperature
- Fever alarm
- 9 memories
- °C / °F switchable
- PLUG & SECURE Probe Cover
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



Patent Nos : ZL 2011 3 0359483.1
ZL 2010 2 0668110.2
M403995

Clinically Accurate

EN 12470-5:2003

This Ear thermometer has been clinically tested in several big teaching hospitals based on the protocol of EN 12470-5:2003 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

Model	Qty per carton	Carton volume
RA600	50 pcs	0.072 cbm/ctn

Storage and Transportation Condition
-25°C~55°C(-13°F~131°F)
≤95% RH, 700~1060 hPa