

Handheld Pulse Oximeter

Model: SA210



- **%SpO2 & Heart Rate**
Backlight LCD display for easy reading %SpO2 & Heart Rate
- **Personal SpO2, Heart Rate, Memory Interval thresholds setting**
- **Shield Design** against ambient light
- **Biocompatibility Protection**
- **Alarm Off**
Device delivers audible warning when measurement result exceeds the SpO2/heart rate thresholds previously set. User can disable the alarm function by pressing the "Alarm off" button.
- **Silent mode**
Disables the beep sound during the measurement by pressing the "mute" button.
- **Up to 288 Memories**
Could be programmed with 5 intervals: 0/5/10/30/60-minute ("0" means no memory setup)
- **Pulse Strength Indicator**



%SpO2 & Heart Rate



High/Low SpO2, HR & Memory Interval Setting



Shield Design against Ambient Light



Biocompatibility Protection



Alarm & Sound Off



Up to 288 Memories



Backlight



Pulse Strength Indication



Protection Pouch

Optional Accessory:



Adult Probe
PA100



Neonatal Probe
PB100



Pediatric Probe
PC100



Disposable Probe
PD100



Tube Probe
PF100

Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

Specification

- Measurement Range: SpO2: 35%~100%; HR: 30~250 bpm
- Precision: SpO2: 70%~100%; $\pm 2\%$; 35%~69% (unspecified); Pulse rate: 30~250bpm; ± 3
- PC-link: Wired
- Memory: Up to 288 Memories
- SpO2 Accessory: Adult Probe Included, Pediatric / Neonatal Probe Optional
- Power Supply: 4" AA" Alkaline

| Model | Qty per carton | Carton volume |
|-------|----------------|---------------|
| SA210 | 10 pcs | 0.039 cbm/ctn |

- Storage temperature: -25°C - 70°C (-13°F - 158°F)
- Relative humidity: 15% - 90% (no condensing)
- Atmospheric Pressure: 700 to 1013 hPa